











Adaptor Plate Cutout Dimensions



ITEM #	OPENINGS	CASE LOT
19190	With two 6 ³ / ₈ " (16.2 cm) inset holes — Fits 78164 (4½ qt/3.9 L)	4
19192	With two 8 ³ / ₈ " (21.3 cm) inset holes — Fits 78184 (7¼ qt/6.9 L)	4
19194	Adaptor ring fits 19192, 19188 and 19198 adaptor plates — Fits 78184 (7¼ qt/6.9 L)	2
72227	With one 8 ³ / ₈ " (21.3 cm) and one 10 ³ / ₈ " (26.4 cm) inset holes — Fits 78184 (7¼ qt/6.9 L) & 78204 (11 qt/10.4 L)	1
19193	With one 4 ⁷ / ₈ " (12.4 cm) and two 6 ³ / ₈ " (16.2 cm) inset holes — Fits 78154 (2½ qt/2.4 L) & two 78164 (4½ qt/3.9 L)	4
72228	With three 6 ¹ / ₂ " (16.5 cm) inset holes — Fits 78164 (4½ qt/3.9 L) — Fits thermoset wells only	1
19188*	With three 8 ³ / ₈ " (21.3 cm) inset holes, fits 78184 (7¼ qt/6.9 L) — For use with items 72050/72051	1
19195	With six 4 ¹ / ₄ " (10.8 cm) holes — Fits 78710 (1¼ qt/1.2 L)	4
19189	With one 10 ¹ / ₂ " (26.8 cm) hole — Fits 78204 (11 qt/10.4 L)	4
19198	With one 8 ³ / ₈ " (21.3 cm) hole — Fits 78184 (7¼ qt/6.9 L)	2

*Fits 4/3 rethermalizers and 4/3 drop-ins only